



CR-39 for Measuring Exposure to Radon (^{222}Rn)



Track Detectors for Radon




Intercast CR-39 detector



Politrack™ instrument


Track Detectors



Physics of track detectors


When an ionizing charged particle passes through a dielectric material the transfer of energy to electrons results in a trail of damaged molecules along the particle's track.

Radiation Detection and Measurements, G. Knoll




The tracks can be made visible with a chemical treatment, called "Etching". The opening of the track is then of about 5-20 μm depending on the type and energy of the hadrons.

B. Dorschel et al. / Radiation Measurements 37 (2003) 563 – 571




Detecting Radon

Alpha particles from ^{222}Rn , ^{214}Po & ^{218}Po




Different shape of tracks left by Radon and Radon daughters (Polonium)


Cross section:



Microscope view:



Measuring Radon



Tracks left by alpha particles in CR-39, from Radon and Polonium, seen with the naked eye.

Measuring Radon





Image of tracks captured, after etching, with the POLITRACK microscope in a CR-39 detector exposed to Radon gas

Track left by alpha particle



Very close look at an alpha particle's track using an Scanning Electron Microscope (SEM)

QUIZ (a short one)



QUIZ (a short one)

- a) Is Radon visible to our naked eye?
- YES of course NO this is impossible
- b) Where could we find Radon in greatest proportions in our house?
- Cellar Terrace
- c) What is the effect of Radon on your body?
- Provokes lung cancer Makes you more beautiful
- d) What is the best technique to avoid accumulation of Radon in your house?
- Ask advice to an expert Ventilate regularly all rooms and cellar
- e) If we cannot see, feel, smell or taste Radon gas, how do we detect it?
- Detect alpha particles from Rn & Po Call the Police

Answers...

- a) Is Radon visible to our naked eye?
- YES of course NO this is impossible
- b) Where could we find Radon in greatest proportions in our house?
- Cellar Terrace
- c) What is the effect of Radon on your body?
- Provokes lung cancer Makes you more beautiful
- d) What is the best technique to avoid accumulation of Radon in your house?
- Ask advice to an expert Ventilate regularly all rooms and cellar
- e) If we cannot see, feel, smell or taste Radon gas, how do we detect it?
- Detect alpha particles from Rn & Po Call the Police

Experiments



*Neil Patrick Harris,
Dr Horrible*

Experiments

- Session 1

- Calibrate CR-39 detectors


- Session 2

- Measure the amount of radon in your house
- Calculate the amount of Radon recorded by the CR-39 detector

Session 1 - Calibration

- a) Using the optical microscope to count the number of track on the CR-39 detectors, left from Radon and Radon daughters for the known Exposure levels.
- b) Calculate the area used for measurements
- c) Calculate the track density (N° of tracks per cm²)
- d) Plot the graph of Track density (N° of tracks per cm²) VS. Radon exposure (Bq*h/m³)
- e) Calculate the **Calibration Factor (C_f)** which is the gradient of the graph (see next slide)

Example of calibration curve




$$\text{Track density} = C_f \times \text{Exposure}$$

Session 2 – Radon measurement

1. Assemble CR-39 detector and holder
2. Expose detectors in a “living space” during period of time T
3. After exposition, detector is etched in conc. NaOH for 60 mins at Mi.am SRL
4. Detectors are read with an optical microscope to count the track density (N° of tracks per cm²)
5. Calculate Exposure using equation 1
6. Calculate Radon Concentration using equation 2
7. CONGRATULATIONS! You are a Radiation Protection Expert!!!


Step 1 – Detector assembly

Take the CR-39 detector



(a)

Put CR-39 in base of holder



(b)


Clip in the cover of the holder



(c)

Special thanks to Barbara Rossi for lending her hands for this presentation

Recap: Where do we find ^{222}Rn ?



- Radon is emanated from the Uranium rich soil or rocks
- Radon can also be found in water
- Radon escapes easily from the ground into air where it decays into its *progeny*
 ^{214}Po , ^{218}Po

World Health Organisation

Step 2 – Radon measurement

This is not my house, these pictures were taken randomly from the internet



Living room



Bedroom



Cellar



Bathroom

The detector should be exposed in a “living space” for several months:

- Bedroom
- Living room
- Study room
- Basement used as playroom
- Cellar (Wine)
- Bathroom

Choose one of those rooms in your house/apartment for the experiment

Step 3 - Etching




Etching is important in order to reveal the tracks formed at the surface of CR-39 after it has been exposed to Radon gas

60 mins in conc. aq. solution of NaOH at 98 °C

Step 4 – Track density (Tracks per cm²)

1. Look at the CR-39 detector under the microscope.
2. Count the total number of tracks you can find on a surface of $1\text{ cm} \times 1\text{ cm} = 1\text{ cm}^2$
3. This number of tracks in 1 cm^2 is the track density
4. Note down this track density for step 5



Step 5 – Calculate Radon Exposure

$$\text{Exposure} = \frac{\rho_{tracks}}{C_f} \quad [\text{Bq} \cdot \text{h} / \text{m}^3]$$

Where ρ_{tracks} is the track density expressed in cm^{-2} and C_f is the calibration factor calculated in Session 1.

Step 6 – Calculate Radon concentration

$$\text{Radon Conc.} = \text{Exposure} / T \quad [\text{Bq}/\text{m}^3]$$

Where *Exposure* is the Radon exposure expressed in Bq*h/m³ and *T* is the exposition time in hours.

NB: Remember to take note of the irradiation time, note the date you start and finish the experiment



JABLONTRON

mi.am



POLITECNICO
DI MILANO

FAU



UNIVERSITY OF
HOUSTON

UOIT

UNIVERSITY OF
WOLLONGONG

Step 7 – Congratulations!!!

You are a
Radiation Protection Expert!!!!

I am joking :)



POLITECNICO
DI MILANO




Applications of CR-39 detector

Widely used for several applications:

- Radon measurement – Environmental dosimetry ★
- Fast neutron dosimetry – Space, Aircrafts and Particle Accelerators
- Thermal neutron dosimetry – Nuclear Reactors
- Cosmic rays detection – Observatories
- High Energy Physics experiments – Nuclear Fusion, Military & Defense

Neutron detection



←→
10 μm

Image of tracks captured, after etching, by POLITRACK in a CR-39 detector (coupled to PMMA) exposed in a quasi-monoenergetic neutron beam.



www.nasa.gov

26/11/2013



www.nasa.gov